

Best practices in the analysis of RNA-seq and ChIP-seq data

27th – 31st, July 2015

University of Cambridge, Cambridge, UK

The quality of a ChIP-seq data

Ines de Santiago

CRUK Cambridge Research Institute

Ines.desantiago@cruk.cam.ac.uk




**UNIVERSITY OF
CAMBRIDGE**




Acknowledgments

- Tom Carroll
 - http://bioconductor.org/help/course-materials/2014/BioC2014/ChIPQC_Presentation.pdf
 - https://github.com/bioinformatics-core-shared-training/ngs-in-bioc/blob/master/Lectures/Lect6b_ChIP-Seq%20Data%20Analysis.pdf
- Shamith Samarajiwa
- Suraj Menon

“Typical” ChIP-Seq Analysis workflow




A good ChIP-seq dataset



What do we want:

- Good quality ChIP-seq enrichment over background

Narrow vs Broad peaks



Nature Reviews | **Genetics**


What can go wrong?

- **The specificity of the antibody**
 - poor reactivity against the intended target
 - cross-reactivity with other DNA-associated proteins.
- **degree of enrichment** achieved in the affinity precipitation step.
- **Biases in the selected fragments:**
 - PCR amplification biases
 - Fragmentation biases

How to quantify ChIP-seq data quality?

- **ChIP-seq guidelines and practices of the ENCODE and modENCODE consortia.**
(Landt et al – *Genome Research* 2012)
- **ChIPQC** – Tom Carroll and Rory Stark (*Diffbind*)
- **ChIPQC** provides workflow to generate metrics per sample/experiment.
- package **SPP** (for UNIX/LINUX)

EVALUATING CHIP-SEQ DATA (QC)



What do we want:

- Good quality ChIP-seq enrichment over background

Outline

- **Distribution of Signal**
 - Visualisation of coverage profiles
 - Fraction of reads in peaks (FRIP)
 - Relative enrichment in genomic intervals (REGI)
 - Signal in blacklists (FRIBL)
 - Dispersion of coverage
- **Clustering of Watson/Crick reads.**
- **Other factors affecting site discovery:**
 - Sequencing depth
 - Duplication rate / library complexity
 - Control sample

Outline

- **Distribution of Signal**



- Visualisation of coverage profiles
- Fraction of reads in peaks (FRIP)
- Relative enrichment in genomic intervals (REGI)
- Signal in blacklists (FRIBL)
- Dispersion of coverage

- **Clustering of Watson/Crick reads.**

- **Other factors affecting site discovery:**

- Sequencing depth
- Duplication rate / library complexity
- Control sample

Visualisation of coverage profiles




Outline

- **Distribution of Signal**
 - Visualisation of coverage profiles
 - Fraction of reads in peaks (FRIP)
 - Relative enrichment in genomic intervals (REGI)
 - Signal in blacklists (FRIBL)
 - Dispersion of coverage
- **Clustering of Watson/Crick reads.**
- **Other factors affecting site discovery:**
 - Sequencing depth
 - Duplication rate / library complexity
 - Control sample

Measuring global ChIP enrichment (FRiP)

- useful and simple first-cut metric for the success of the immunoprecipitation
- Good quality TF > 5%

Example output from ChIPQ package:




Outline

- **Distribution of Signal**
 - Visualisation of coverage profiles
 - Fraction of reads in peaks (FRIP)
 - Relative enrichment in genomic intervals (REGI)
 - Signal in blacklists (FRIBL)
 - Dispersion of coverage
- **Clustering of Watson/Crick reads.**
- **Other factors affecting site discovery:**
 - Sequencing depth
 - Duplication rate / library complexity
 - Control sample

Enrichment in genomic intervals

- Plot relative enrichment of reads in annotated regions.

Example output from ChIPQ package:




Outline

- **Distribution of Signal**
 - Visualisation of coverage profiles
 - Fraction of reads in peaks (FRIP)
 - Relative enrichment in genomic intervals (REGI)
 - Signal in blacklists (FRIBL)
 - Dispersion of coverage
- **Clustering of Watson/Crick reads.**
- **Other factors affecting site discovery:**
 - Sequencing depth
 - Duplication rate / library complexity
 - Control sample

Signal in blacklists

- Encode empirically identified regions that showed anomalous and near-universal artefact signal
- Represent around 0.5% of genome.
 - Various reasons e.g. chromatin accessibility, repeats
- Can account for high proportion of total signal (> 10%)
 - Enriched for duplicate and multi-mapping reads
- Adversely affect fragment length calculations and in thus any analyses that require these e.g. peak calling

Example output from ChIPQ package:




Outline

- **Distribution of Signal**
 - Visualisation of coverage profiles
 - Fraction of reads in peaks (FRIP)
 - Relative enrichment in genomic intervals (REGI)
 - Signal in blacklists (FRIBL)
 - Dispersion of coverage
- **Clustering of Watson/Crick reads.**
- **Other factors affecting site discovery:**
 - Sequencing depth
 - Duplication rate / library complexity
 - Control sample

Dispersion of coverage

- Depth of signal: number of fragments at a genomic location.
- Expectation is that for an enriched ChIP sample, depth should show inequality in dispersion across the genome
- Build global profile of signal depth
 - Measure number of base pairs with given depth of signals.
 - Normalise to total number of reads to compare samples



Depth	Base Pairs
1	3
2	4
3	3
5	3
6	4
7	3
8	26

Dispersion of coverage

- Global signal profile “histogram”
- Enriched (ChIP) libraries show higher number of bases at greater depths.
- Profile for inputs (no enrichment) drops off more quickly
- Gap between sample and input indicates enrichment

Example output from ChIPQ package:



Metric for dispersion of coverage: SSD


- SSD: Standardised Standard Deviation of coverage

- Metric for assessment of dispersion coverage developed in htseqtools package

$$SSD = \frac{SD}{\sqrt{n}}$$

- Provides measure of pile-up across genome
 - High for samples with enriched regions (ChIP)
 - Low for samples with uniform coverage (input)

Example output from ChIPQ package:





Outline

- **Distribution of Signal**
 - Visualisation of coverage profiles
 - Fraction of reads in peaks (FRIP)
 - Relative enrichment in genomic intervals (REGI)
 - Signal in blacklists (FRIBL)
 - Dispersion of coverage
- **Clustering of Watson/Crick reads.**
- **Other factors affecting site discovery:**
 - Sequencing depth
 - Duplication rate / library complexity
 - Control sample

Clustering of Watson/Crick reads

How to make a cross-correlation plot:




<http://www.nature.com/nmeth/journal/v6/n4/images/nmeth.f.247-F2.jpg>
<http://www.bloodjournal.org/content/124/25/3719>

Clustering of Watson/Crick reads


- Fragment length can be estimated from data:
 - **Cross-correlations** - Correlation of reads on positive and negative strand after successive read shifts
 - **Cross-coverage** - Coverage of reads on both strand after successive shifts of reads on one strand. Total area covered by reads will be reduced after shifting
- These provide useful QC metrics

How to make a cross-correlation plot:



Clustering of Watson/Crick reads

- Cross-correlation/Cross-coverage score plots are useful for QC
- Blacklisted regions strongly contribute to read length cross-coverage peak
 - Small to non-existent peaks are seen in failed ChIPs and inputs
- ChIPQC metrics:
 - $FragCC = CC_{fragmentlength}$
 - $RelCC = FragCC / CC_{readlength}$




Outline

- **Distribution of Signal**
 - Visualisation of coverage profiles
 - Fraction of reads in peaks (FRIP)
 - Relative enrichment in genomic intervals (REGI)
 - Signal in blacklists (FRIBL)
 - Dispersion of coverage
- **Clustering of Watson/Crick reads.**
- **Other factors affecting site discovery:**
 - Sequencing depth
 - Duplication rate / library complexity
 - Control sample

Sequencing Depth

Peak counts depend on sequencing depth.



Sequencing Depth: guidelines

Sharp peaks (TFs)

10M reads

2M worms and flies

Broad Peaks (Histones)

20M reads mammalian genomes

5M worms and flies

Outline


- **Distribution of Signal**
 - Visualisation of coverage profiles
 - Fraction of reads in peaks (FRIP)
 - Relative enrichment in genomic intervals (REGI)
 - Signal in blacklists (FRIBL)
 - Dispersion of coverage
- **Clustering of Watson/Crick reads.**
- **Other factors affecting site discovery:**
 - Sequencing depth
 - Duplication rate / library complexity
 - Control sample

Library complexity (Duplicates)

A



Typical ChIP-seq peak



Low-complexity ChIP-seq peak

Library complexity (Duplicates)

- Duplicates can be artefacts
- PCR bias: certain genomic regions are preferentially amplified
- Low initial starting material
 - Overamplification -> artificially enriched regions
 - Compounded by PCR bias
- Duplicates can also be 'legitimate'
 - In highly efficient enrichments
 - In deeply sequenced ChIPs (Duplication rate increases with sequencing depth)
- Removing these duplicates limits the dynamic range of ChIP signal
 - Max signal for a base is $(2 * \text{read length}) - 1$

Library complexity (Duplicates)

- So what to do about duplicates?
- Keep in mind enrichment efficiency and read depth
- Thumb-rules
 - Remove duplicates prior to peak calling (some peak callers do this by default)
 - Keep duplicates for differential binding analysis
- A more objective approach:
 - htSeqTools package
 - Estimate duplicate numbers expected for sequencing depth using negative binomial model and attempt to identify significantly anomalous duplicate numbers.

Library complexity (Duplicates)

- Duplication rates are a useful QC metric
 - $(\text{Duplicate reads} / \text{Total Mapped Reads}) * 100$
 - Expected to be low ($< \sim 1\%$) for inputs
- Non-Redundant Fraction (NRF)
 - $\text{Unique Reads} / \text{Total Mapped Reads}$
 - ENCODE guidelines:
 - NRF ≥ 0.8 for 10M reads

Outline


- **Distribution of Signal**
 - Visualisation of coverage profiles
 - Fraction of reads in peaks (FRIP)
 - Relative enrichment in genomic intervals (REGI)
 - Signal in blacklists (FRIBL)
 - Dispersion of coverage
- **Clustering of Watson/Crick reads.**
- **Other factors affecting site discovery:**
 - Sequencing depth
 - Duplication rate / library complexity
 - Control sample

Control sample

- Use of controls **highly** recommended
- **Input DNA**
 - popularly used
 - controls for CNVs, sequencing biases, fragmentation and shearing biases
- **IgG**
 - as with input but also controls for non-specific binding
 - but introduces new biases
- **Controls required for**
 - different types of samples (e.g. Cell lines, mice, patients)
 - different treatment groups / experimental conditions

Control sample

A Normal input



B Input with strong “Sono-Seq” effect




Fig. S1

PEAK CALLING


Peak Calling:

Which Peak Caller to Use?

- Transcription factor peaks: **MACS** is very popular
- For histone marks with spanning longer regions, **Sicer** is recommended
 - MACS can be used by tweaking parameters
- Several peak callers in R/Bioconductor
 - e.g SPP, TPIC, BayesPeak
 - Not really considered gold-standard (other than SPP)
 - Often impractical: memory hungry and slow

Peak Calling: MACS

- Sliding window run across genome
- Peak height in window compared to that in windows of surrounding regions in control





- Statistical significance of peak estimated by using Poisson distribution
 - $-\log_{10}(\text{pvalue})$ reported as peak score
- FDR calculated by calling peaks in control over sample

Peak Calling: Post-peak QC

Example output from ChIPQ package:

- Peak profile plots
 - Mean read density at positions relative to peak summits
 - Input profiles should be flat
- Fraction of Reads in Peaks (FRIP)
 - Reads in peaks/Total mapped reads
 - Analogous to signal to noise ratio



ChIP-Seq Practical

Working with ChIP-Seq Data in R/Bioconductor

chipqc_sweave.pdf